

Otto Bock®

Avantgarde T



The „T“ in Avantgarde T stands for top-of-the-line, a wheelchair for active use, packed with functions. Attractively designed to master demanding everyday situations. It adapts to your needs, for example, with the new ergonomic back angle adjustment of -6° to +15°. The height of the pivot point is also adjustable and can be adapted individually to seat cushions of any thickness. The new drive wheel adapter provides 13 positions for the rear seat height, from 37.5 to 53 cm, as well as 6 adjustment possibilities for variously positioned centres of gravity. Its design is a match for the functional demands placed on it: from new side panels to innovative caster wheel forks and optimised adjustment belts for the back as standard, making the Avantgarde T as a whole very convincing.

Form follows function – It's T time!



Avantgarde T in Detail

Anatomically adjustable backrest with new belt system • Height-adjustable backrest standard • New aluminium side panels • New footrest adjustable for height and angle • Precisely adjustable drive wheel position thanks to optimised adapters • Angle of bend of the front frame is selectable • Extensive options, including new optional ergonomic backrest with a physiologically adjustable pivot and angle adjustment from -9° to $+15^{\circ}$.

Technical Data

Seat width	32-50 cm
Seat depth	36-48 cm
Back height	25-50 cm
Front seat height	40-56 cm
Rear seat height	37-53 cm
Max. user weight	125 kg
Weight	from 11.5 kg
Seat tilt	adjustable up to 15°

Highly rigid aluminium frame with double cross braces

New aluminium footrest with optimised locking function and functional design